


## Sample ID panel workflow

This protocol contains a multiplex PCR step for the simultaneous production of 50 amplicon targets in a single tube and an indexing step for the addition of dual indexed adapters,

### STEP 1

### STEP 2



## Required components

|                           |                         |
|---------------------------|-------------------------|
| Sample ID panel plex (5X) | included in the package |
| pxlence indexes (17X)     |                         |
| pxlence PCR mastermix     |                         |

## Preparations

Create the following cycling conditions in your preferred thermal cycler.

### Multiplex PCR

|        | Temperature (°C) | Time (s) | Number of repetitions |
|--------|------------------|----------|-----------------------|
| Step 1 | 98               | 180      |                       |
| Step 2 | 98               | 15       | repeat 40 times       |
| Step 3 | 60               | 240      |                       |
| Step 4 | 10               | ∞        |                       |

### Indexing PCR

|        | Temperature (°C) | Time (s) | Number of repetitions |
|--------|------------------|----------|-----------------------|
| Step 1 | 98               | 15       |                       |
| Step 2 | 60               | 240      | repeat 5 times        |
| Step 3 | 10               | ∞        |                       |

## Protocol

### Multiplex PCR

1. Create the multiplex PCR reaction by mixing the components as follows in a 0.2 ml eppendorf :

|                                  | Volume (µL) |    |
|----------------------------------|-------------|----|
| SsoAdvanced Preamp Supermixv(2X) | 7.5         |    |
| Sample ID panel plex (5X)        | 3           |    |
| Template DNA (1-5 ng)            | 1 to 4.5    |    |
| Water                            | 0 to 3.5    |    |
| Total reaction volume :          |             | 15 |

2. Place the eppendorf in a thermocycler and run the **Multiplex PCR** thermal cycling conditions.

### Indexing PCR

3. Take the eppendorf out of the thermal cyclcr when the **Multiplex PCR** thermal cycling conditions have been completed.
4. Add 1 µL of the selected 17X i7-index solution and 1 µL of the selected 17X i5-index solution to the eppendorf.
5. Vortex carefully or pipet 10 times up and down.
6. Place the eppendorf in a thermocycler and run the **Index PCR** thermal cycling conditions.
7. When the **Index PCR** thermal cycling conditions have been completed, take the eppendorf out of the thermal cyclcr and store it at -20°C until sequencing.