

## Rainbow probes: universal detection of DNA for digital PCR

### What are Rainbow™ probes?


Rainbow probes are patent-pending fluorescent nucleic acid detection probes that eliminate the need for designing and optimizing target-specific probes. Here's how they work:

- **Modify existing primers**

We transform one of your existing PCR primers (from a pair) into a Rainbow™ primer by modifying the oligonucleotide. This modified primer allows the Rainbow™ probe to detect the amplicon.

- **Multiplex capability**

With 2 to 8 Rainbow™ probes (depending on your digital PCR instrument), you can quantify multiple targets simultaneously (multiplex reaction). Perfect for applications requiring numerous probes or small probe quantities.



### Benefits

- + **Time saving**

No design and optimization of probe sequences required; fast as probes are delivered off the shelf.

- + **Cost efficient**

Use the same probe for different targets, minimizing unused stock.

- + **High sensitivity (1)**

Rainbow probes are 60X concentrated, enabling efficient and sensitive multiplexing by accommodating large sample input volume in the reaction.

- + **High sensitivity (2)**

No sequence space needed for a probe, allowing ultrashort amplicons, increasing analytical sensitivity when working with fragmented nucleic acids from e.g. FFPE tissue or liquid biopsy.


## Available Rainbow™ probes

probe name	primer tag	fluorophore	emission max (nm)	compatible instruments
RBO-01-A-425	A	ATTO 425	484	Digital LightCycler, DQ24
RBO-01-B-FAM	B	FAM	517	QX200, QX ONE, QX600, QIAcuity, naica 3, naica 6, Nio, Digital LightCycler, DQ24, QuantStudio Absolute Q
RBO-01-C-YAK	C	Yakima Yellow	550	naica 6, Nio
RBO-01-C-HEX	C	HEX	559	QX200, QX ONE, QX600, QIAcuity, naica 6, Digital LightCycler, DQ24, QuantStudio Absolute Q
RBO-01-A-550	A	ATTO 550	574	QIAcuity, naica 6, Nio, QuantStudio Absolute Q
RBO-01-D-TEX	D	Texas Red	603	QIAcuity, Digital LightCycler
RBO-01-D-ROX	D	ROX	604	QX600, QIAcuity, naica 6, Nio, DQ24
RBO-01-A-590	A	ATTO 590	621	QX600
RBO-01-F-CY5	F	Cy5	670	QX ONE, QX600, QIAcuity, naica 3, naica 6, Nio, Digital LightCycler, DQ24, QuantStudio Absolute Q
RBO-01-G-C55	G	Cy5.5	703	QX ONE, QX600, Digital LightCycler, DQ24
RBO-01-G-700	G	ATTO 700	715	naica 6, Nio


## Representative 1-color plots

Below are representative 1-color singleplex results using different Rainbow™ probes on different digital PCR instruments. Each time, a BRAF V600E wild type assay was used with sheared human genomic DNA as input (0.7 copies per partition ( $\lambda$ ), resulting in an equal number of positive and negative partitions).


Roche Digital LightCycler



Bio-Rad QX200



Qiagen QIAcuity



## Contact

URL : [www.pxlence.com](http://www.pxlence.com)  
 email : [info@pxlence.com](mailto:info@pxlence.com)

LinkedIn : [www.linkedin.com/company/pxlence](http://www.linkedin.com/company/pxlence)

[more info:](#)

