

## PCR based target enrichment for variant confirmation, gene panels and multiplex PCR sample tracking in a whole exome sequencing workflow


### Product portfolio

Our offering comprises **best-in-class PCR assays** for targeted resequencing of all human canonical exons of protein coding genes (both Sanger and massively parallel sequencing)


#### WGS/WES variant validation

 **1 million PCR assays**, of which over 10,000 have been validated for variant confirmation

#### Sample ID panel

 An **optimized 50-SNP multiplex panel** to identify your samples in a single PCR reaction

#### NGS gap filling

 Extend your preferred enrichment method with our assays to **fill low coverage regions**


All our assays have been thoroughly validated *in silico*, resulting in an off-the-shelf **PCR success rate of >97%**. The optimized design parameters allow **uniform PCR conditions** while ensuring **uniform sequencing coverage**.

- ✓ Avoid SNPs in primer annealing sites
- ✓ Minimize secondary structures
- ✓ Maximize assay specificity
- ✓ Optimize GC content, annealing temperature and other parameters

### Single-tube 2-step multiplex PCR sample ID sequencing panel

#### Characteristics

- ⊙ Easy-to-use protocol (automatable)
- ⊙ **1 h hands-on time**, 5 h total time
- ⊙ 50 polymorphic SNPs, including 5 gender markers
- ⊙ **No clean-up** prior to sequencing
- ⊙ **Only 4 pipeting steps**
- ⊙ Works with high-quality, FFPE and liquid biopsy samples
- ⊙ Compatible with Illumina sequencing platforms




### Best-in-class assays for WGS/WES variant validation

Almost **1 million PCR assays** for targeted resequencing of all human canonical exons of protein coding genes

Our two available assay databases - a long (350-750 bp) and short (125-275 bp) assay-set - cover over **97.99%** and **98.71%** of the **human canonical exome**, respectively. In approximately **94.01%** of these assays, **no SNPs** are present in the primer annealing sites. For the remaining assays, SNP(s) are located outside the critical 5 bp 3' region. Over **77%** of the assays, **no aspecific annealing sites** harboring less than 3 mismatches can be identified, resulting in a minimal non-specific product generation.

Our short assay delivery times allows **WGS/WES variant identification to validation** in approximately **one week**.



### Contact

www.pxlence.com  
 frauke.coppieters@pxlence.com