

## PCR-grade water

Parameter	Specification
<i>Measured on raw material before filtration</i>	
Bioburden	Conform Ph. Eur. Monograph 0008
Conductivity	Conform Ph. Eur. Monograph 0008
TOC	<0.5 ppm
Nitrates	<0.2 ppm NO <sub>3</sub>
Appearance	Clear, colorless fluid


### *Measured on final product*

Sterility test (performed according to the current Ph. Eur. 2.6.1 / USP)	No growth
RNase activity	None detected
DNase activity	None detected
Endotoxins	<0.25 EU/mL

## Screw-capped polypropylene tubes

Volume: 1.5 mL

Quantity: 72



## PETG bottles

Volume: 100 mL

Quantity: 1



Volume: 1000 mL

Quantity: 1

